

# CERTIFICATE OF ANALYSIS

Document Code: QO1N4501-R12      Implementation Date: 2023.09.01      Version: T/00      Page: 1/1

## 1 Product Information

Product Name	Rapid SARS-CoV-2/Flu A+B/RSV Antigen Test Card		
Lot Number	24071501	Mfg. Date	2024.07.15
Quantity	195020 Tests	Exp. Date	2026-01-14

## 2 Report Contents

Inspection Item	Standards	Result
Physical Inspection	Components should be complete and intact. Packaging should be undamaged. Liquid should be clear with no leakage. Labels should be firmly attached with clear and correct contents.	Pass
Limit of detection	Test results for LoD control S1 (SARS-CoV-2) should be positive for SARS-CoV-2 and negative for FluA, FluB and RSV.	Pass
	Test results for LoD control S2 (FluA) should be positive for FluA and negative for SARS-CoV-2, FluB and RSV.	Pass
	Test results for LoD control S3 (FluB) should be positive for FluB and negative for SARS-CoV-2, FluA and RSV.	Pass
	Test results for LoD control S4 (RSV) should be positive for RSV and negative for SARS-CoV-2, FluA and FluB.	Pass
Negative agreement	All test results for negative controls should be negative for SARS-CoV-2, Flu A, Flu B and RSV.	Pass
Positive agreement	Test results for positive controls P1 and P2 should be positive for SARS-CoV-2 and negative for FluA, FluB and RSV.	Pass
	Test results for positive controls P3 and P4 should be positive for FluA and negative for SARS-CoV-2, FluB and RSV.	Pass
	Test results for positive controls P5 and P6 should be positive for FluB and negative for SARS-CoV-2, FluA and RSV.	Pass
	Test results for positive controls P7 and P8 should be positive for RSV and negative for SARS-CoV-2, FluA and FluB.	Pass

The lot (check one)      ☒ passes    ☐ fails      QC test

Authorized signature: \_\_\_\_\_

Date: 2024.07.18